



CO2 Euthanasia programmer


This programmer is dedicated to labs rodents euthanasia and allows a CO2 induction session following a strict and secured protocol.

Product features

A progressive anesthesia period, saturation with CO2, dwell period and extraction of CO2. This device is an electronic controlled pressure circuitry for a better comfort of animals.

Technical specifications

Only dedicated for rodents (except new born)
Operator safety
Easy cleaning
Automatic protocols
Lime CO2 filter
Optimized CO2 consumption
Electronic pressure sensor



Compliant to European directive
2010/63/UE