



CO2 Euthanasia Housing

This housing is dedicated to labs rodents euthanasia and allows a CO2 induction session following a strict and secured protocol.

Product features

A progressive anesthesia period, saturation with CO2, dwell period and extraction of CO2. This device is an electronic controlled pressure circuitry for a better comfort of animals.

Technical specifications

- Only dedicated for rodents (except new born)
- Full stainless mobile housing
- Operator safety
- Easy cleaning
- Automatic protocols
- Electric opening safety
- Optimized CO2 consumption
- Electronic pressure sensor

