



GenieTouch™

Infusion / withdrawal dual syringe pumps

- Multi-syringe capability
- Push/Pull capability
- Touch screen user interface
- Automatic positioning/stalling detection

■ Syringe specifications

From 0,5 µl to 100 ml.
Up to 4 simultaneous syringes
Empty / full point setting.

■ Injection modes

Constant rate, time dispensing on volume dispensing, ramps, pulsed, stepped.
From some nanoliter/min up to some ml/min.

■ Advanced features

Programmable alarms
Pusher block stall detection
Automatic calibration
Autoreverse mode
Adjustable linear force
External trigger
USB interface

