



MiniHub

Ultra-compact gaseous anesthesia station

- For small animals
- Low isoflurane consumption
- Designed for all type of hubs
- Flow regulation for each output
- Easy configurable touch screen
- Indication of filter saturation



■ Dimensions

Depth : 41 cm - Width: 50 cm- Height: 33 cm

■ Security

Automatic detection of opening/closing of the box
Automatic detection of animal presence on the mask
Electronic valves and flowmeters
1 general control panel with touch screen
1 flowmeter O² or AIR
Induction chamber temperature regulation

■ Peripherals

1 ISOTEC3 or ISOTEC7 vaporizer with secured filling system
1 or 2 induction chamber
1 or 2 automatic mask for rat or mice
USB port – protocols memorization

■ Scavenging system

- Regulated scavenger for mask
- Double safety system : internal and external scavenging
- Easy removable carbon filter rack
- Visual alarm for filter saturation

