



MiniTag

Gaseous anesthesia station MiniTAG

- For small animals
- Fully integrated and automatic device
- 2 easy configurable workstations
- 2 carrier gas mixing system
- Proactive embedded scavenging system
- One Touch interface

■ Dimensions

Depth: 80 cm.

Usable workstation depth: 70 cm.

Width: 66cm

Height: 90 cm.

■ Mecanics

4 directional wheels

Stainless liquid rack with purge tap for easy cleaning.

■ Gas circuitry

Embedded electronic valves and flowmeter

Air compressor embedded

1 touchscreen integrated command panel

1 temperature sensor

2 flowmeters O² & AIR

■ Peripherals

1 ISOTECH7 vaporizer with securized filling system

1 induction chamber

1 mask for rat or mice

2 perforated stainless plates with 1 heating plate.

■ Scavenging system

- Easy removable carbon filter rack

- Visual alarm for filter saturation

- Embedded dust filter.

