



NUVO

Oxygen concentrator

Product features

This concentrator has been validated for use in labs research centers. It can be connected to any anesthesia station.

Features

- Small size device
- Easy access to filters
- Flow from 0,125 to 5L/min
- Mobile device
- Low cost use
- Very silent device <40 dB

Technical specifications

- Power supply 230V, 50 Hz, 290W
- Storage temperature : - 20 à + 60°C
- Using temperature : + 5 à + 40°C
- Weight : 13,6 Kg
- Dimensions : 35 x 58 x 22 cm [l x h x p]
- Norms : ISO 8359, EN 60601-1
- Class CE : II b CE 0413
- Concentration O2 : at 2L/min : >90%,
at 5L/min : 90%
[+6,5% : +3%]



Not compatible with MRI environment