



Tabletop Anesthesia Station

Dedicated for labs with Air and/or Oxygen network, this Tabletop is very easy to use.

Delivered with induction chamber and masks, it is mainly used for rodents, but can also be used for birds and all others small animals

■ Product Features

- Can support many anesthesia interfaces: for hub, MRI, scanner, bioluminescent imager, surgery, etc,
- Can support many add-ons: stereotaxic mask, warwing systems,...
- Pneumatic scavenger for halogen protection

■ Technical specifications

Gas alimentation

- Stainless tabletop structure.
- Oxygen flowmeter (0-14 L/min).
- Air flowmeter (0-14 L/min).
- Overpressure safety valve.
- TEC3 vaporizer keyFill™ secured filling system. Can support : Isoflurane, Sevoflurane.
- Distribution block for up to 6 gas connections.

Induction

- one induction chamber: standard or autoclavable
- Security system to switch on the scavenger when chamber opening.

Anesthesia interfaces

- Bain circuitry
- Mask depending of your animal race (rat / mice / rabbit / other).
- Experimentation Support for surgery, heating system, etc.

Scavenging system

- Pneumatic scavenger with two independent channels. One channel for induction chambers (with pressure regulation, one channel for masks (with flow meter regulation).
- F AIR canisters to capture halogens.

■ Accessories / options

- Other vaporizers available: TEC7, second-hand TEC3.
- Induction chambers available : standard or stainless.
- Temperature regulation systems (water or electronics)
- Possible flow meters N2O, CO2.

