



Tag



Gaseous anesthesia station

- For small animals
- Fully integrated and automatic device
- 3 easy configurable workstations
- Up to 4 carrier gas mixing system
- Proactive embedded scavenging system
- Touch screen interface

■ Dimensions

Depth: 80 cm.

Usable workstation depth : 70 cm.

Width: 120cm upgradable to 140 or more

Height: 90 cm.

■ Mecanics

4 directional wheels

Stainless liquid rack with purge tap for easy cleaning.

■ Gas circuitry

Embedded electronic valves and flowmeter

Air compressor embedded

4 touchscreen integrated command panel

1 temperature sensor

Up to 4 flowmeters O² & AIR (and CO₂ – N₂O)

■ Peripherals

1 ISOTECH7 vaporizer with securized filling system

1 induction chamber

1 mask for rat or mice

2 perforated stainless plates with 1 heating plate.

■ Scavenging system

- Easy removable carbon filter rack

- Visual alarm for filter saturation

- Embedded dust filter.