



## Teg

### Euthanasia CO2 mobile station

- Only for rodents (except new born)
- Fully automatic
- 2 carrier gas mixing system
- Easy touch screen interface

#### ■ Dimensions

Depth : 80 cm.

Usable workstation depth : 70 cm.

Width : 66cm.

Height : 90 cm.

Weight : 70kg.

#### ■ Mecanics

4 wheels table with 2 brakes.

Stainless housing can receive up to 1 type 4S cage.

Red translucent PMMA Top.

TRESPA workstation.

#### ■ Gas circuitry

Electronic valves and flowmeters

1 embedded air

Compressor 46 dB for 50L/min

Programmable cycles

#### ■ Peripherals

CO2 bottle support water recuperation system

Lime filter in a rack for 5 liters

Low consumption : 1 rack for 12 cycles

