



Therm250

Body temperature regulation device for small animals

- Rectal probes
- Burning safety
- ULTRAFRESH® pad treatment



COMPATIBLE
MRI



■ Features

To regulation from 30° to 45°C.
Audible burning
Safety alarm
Temperature memorization

■ Technical specifications

Double temperature safety : one rectal probe + one pad probe
to protect the animal.
Pad with special treatment : anti-microbial agent
anti-fungal agent

■ Options

- THERM250
Delivered with 250 x 150 mm heating pad
- THERMTAB
Delivered with surgical heating plate 250 x 100mm
- THERM500
Optional 500 x 500 mm heating pad